

## ARTICLE

### Eyewash bottle Barikos - 620 ml - excluding bracket

Number: 4.25.0040

## SPECIFICATIONS

Measurements: 258 (l) x 92 (w) x 55 (h) mm

Weight: 1000,00 kg

Made in: Germany

HS Code: 39233010



## DESCRIPTION

The Barikos eyewash bottle KS contains pH-neutral water-based rinsing fluid with which you can rinse the eye when it has come into contact with harmful substances. This makes this eye wash an indispensable product in places where chemicals or mechanical equipment are used.

The Barikos eyewash bottle contains 620 ml of sustainable water without chemical additives and can be used once. In environments where there is a high risk of harmful substances getting into your eye, this eye wash must be accessible within 6 minutes. It is recommended to rinse for 15 minutes.

The advantage of the Barikos eyewash bottle is that it can be used both vertically and horizontally, so the liquid is released when the bottle is turned upside down and placed on the eye or when it is squeezed when it is held horizontally against the eye.

Instructions for use in case of emergency:

Remove the protective cap from the Barikos eyewash bottle;  
Bend your head over the bottle or hold the eyewash bottle above your eye;

Keep the eye open with your thumb and index finger;  
Press the spray button gently against the eye socket and rinse the opened eye by repeatedly squeezing the eyewash bottle.

Always consult a doctor after eyewash.

The eyewash bottle can be used yourself, but in the event of a serious eye accident it may be necessary to call in someone else for first aid and to assist with the eyewash.

## SPARE PARTS / COMPLEMENTARY GOODS

4.25.0041 Bracket Steel - Eyewash Bottle Barikos

4.25.0043 Bracket / wall holder with dust-proof housing - Eyewash Bottle Barikos