

ARTICLE

EEBD Semmco Marine (15min) - MED/SOLAS

Number: 4.10.0520

SPECIFICATIONS

Measurements: 210 (l) x 100 (w) x 230 (h) mm

Weight: 2,00 kg

Made in: United Kingdom

HS Code: 90200000



DESCRIPTION

EEBD Semmco Marine (15min) - MED/SOLAS

Semmco Marine 15 EEBD – 15 minute duration for military, maritime and offshore applications. The Semmco EEBD uses potassium superoxide (KO₂) to provide oxygen to the user. As the EEBD is pulled from its carry case, a chlorate candle provides oxygen for the user's first breath. As the user exhales into the bag, moisture from their exhaled breath activates the KO₂ that then 'scrubs' the CO₂ and liberates oxygen. Oxygen continues to be generated while the EEBD is being used.

The Semmco EEBD is contained in a rigid case that is discarded once in use and the Semmco EEBD Marine fits into existing storage racks. The case features both an exterior indicator showing the validity of the set at a quick glance, as well as pictorial instructions for use. For industrial applications the set can be belt mounted using the belt clip accessory.

The Semmco EEBD should be removed once a safe area has been reached.

Features & Benefits

- Quick and easy to don – comfortable to wear
 - Automatic quick start mechanism - instant O₂ supply
 - Lightweight (1.3kg)
 - No face fit test required
 - Easily adjustable neck strap
 - Clear face shield for panoramic vision
 - High visibility breathing bag
 - Ideal for confined spaces - no bulky cylinder to carry
 - Maintenance free for up to 15 years (visual inspections recommended)
 - Clear visual indicator on case to check viability
 - Robust, environmentally sealed product packaging featuring side mounted shoulder strap
 - Fits existing EEBD stowage racks on board ships
 - Belt mounted waist clip for industrial applications
- 3072

CERTIFICATIONS

- MED
- SOLAS
- UKCA Approval
- CE